



Product Datasheet

BCL2L10 antibody (orb686392)


Description	BCL2L10 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Bcl-2-like protein 10 protein (19-127AA)
Target	BCL2L10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HD36
Dilution Range	IHC-P: 1: 20-1: 200, IF: 1: 20-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human li...