

Product Datasheet

Ace antibody (HRP) (orb686433)

Description	Ace antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Rat Angiotensin-converting enzyme protein (1046-1263AA)
Target	Ace
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P47820
Expiration Date	12 months from date of receipt.