



Product Datasheet

EIF2AK3 antibody (orb686440)


Description	EIF2AK3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Eukaryotic translation initiation factor 2-alpha kinase 3 protein (115-330AA)
Target	EIF2AK3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZJ5
Dilution Range	IHC-P: 1: 200-1: 500, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sm...



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human ov...