

## Product Datasheet

### NEURL1B antibody (FITC) (orb686466)

|                        |                                                                            |
|------------------------|----------------------------------------------------------------------------|
| <b>Description</b>     | NEURL1B antibody (FITC)                                                    |
| <b>Species/Host</b>    | Rabbit                                                                     |
| <b>Reactivity</b>      | Human                                                                      |
| <b>Conjugation</b>     | FITC                                                                       |
| <b>Immunogen</b>       | Recombinant Human E3 ubiquitin-protein<br>ligase NEURL1B protein (1-275AA) |
| <b>Target</b>          | NEURL1B                                                                    |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol             |
| <b>Form/Appearance</b> | Liquid                                                                     |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid<br>repeated freeze.           |
| <b>Note</b>            | For research use only                                                      |
| <b>Isotype</b>         | IgG                                                                        |
| <b>Clonality</b>       | Polyclonal                                                                 |
| <b>Uniprot ID</b>      | <b>A8MQ27</b>                                                              |
| <b>Expiration Date</b> | 12 months from date of receipt.                                            |