

Product Datasheet

RUBCN antibody (FITC) (orb686562)

Description	RUBCN antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Run domain Beclin-1-interacting and cysteine-rich domain-containing protein (518-733AA)
Target	RUBCN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92622
Expiration Date	12 months from date of receipt.