

Product Datasheet

WDR35 antibody (Biotin) (orb686591)

Description	WDR35 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human WD repeat-containing protein 35 protein (954-1181AA)
Target	WDR35
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P2L0
Expiration Date	12 months from date of receipt.