

Product Datasheet

Collagen VI antibody (Biotin) (orb686639)

Description	Collagen VI antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Collagen alpha-1(VI) chain protein (64-80AA)
Target	COL6A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12109
Expiration Date	12 months from date of receipt.