


Product Datasheet

ITGA8 antibody (orb686688)

Description	ITGA8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence from Human Integrin alpha-8 protein (114-131AA)
Target	ITGA8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53708
Dilution Range	IHC-P: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...