



Product Datasheet

PDE1B antibody (orb686804)

Description	PDE1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1B protein (1-277AA)
Target	PDE1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01064
Dilution Range	WB: 1: 500-1: 2000, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells using PDE1B antibody.



Western blot analysis of MCF-7 whole cell lysate using PDE1B antibody.