


Product Datasheet

NKAPL antibody (orb686868)

Description	NKAPL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence from Human NKAP-like protein (114-123AA)
Target	NKAPL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5M9Q1
Dilution Range	WB: 1: 500-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...