

Product Datasheet

CAPRIN2 antibody (FITC) (orb686890)

Description	CAPRIN2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Caprin-2 protein (83-100AA)
Target	CAPRIN2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IMN6
Expiration Date	12 months from date of receipt.