


Product Datasheet

TNFAIP8L2 antibody (orb686908)

Description	TNFAIP8L2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 2 protein (1-184AA)
Target	TNFAIP8L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P589
Dilution Range	IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HepG2
cel...