

Product Datasheet

TNFAIP8L2 antibody (HRP) (orb686909)

Description	TNFAIP8L2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 2 protein (1-184AA)
Target	TNFAIP8L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P589
Expiration Date	12 months from date of receipt.