


Product Datasheet

CTHRC1 antibody (orb686928)

Description	CTHRC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Collagen triple helix repeat-containing protein 1 protein (31-243AA)
Target	CTHRC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96CG8
Dilution Range	IHC-P: 1: 200-1: 500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...