


Product Datasheet

TCF7L1 antibody (orb687012)

Description	TCF7L1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence from Human Transcription factor 7-like 1 protein (80-99AA)
Target	TCF7L1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCS4
Dilution Range	WB: 1: 500-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...