

Product Datasheet

CACNA1B antibody (FITC) (orb687038)

Description	CACNA1B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Voltage-dependent N-type calcium channel subunit alpha-1B protein (2132-2339AA)
Target	CACNA1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00975
Expiration Date	12 months from date of receipt.