

Product Datasheet

ZFYVE9 antibody (HRP) (orb687161)

Description	ZFYVE9 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Zinc finger FYVE domain-containing protein 9 protein (683-926AA)
Target	ZFYVE9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95405
Expiration Date	12 months from date of receipt.