

Product Datasheet

VPS36 antibody (FITC) (orb687210)

Description	VPS36 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Vacuolar protein-sorting-associated protein 36 protein (87-105AA)
Target	VPS36
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VN1
Expiration Date	12 months from date of receipt.