



Product Datasheet

APCDD1 antibody (orb687212)


Description	APCDD1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protein APCDD1 protein (22-288AA)
Target	APCDD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8J025
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 20-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human he...



Immunohistochemical staining of human pa...



Western blot analysis of 293T whole cell...