

Product Datasheet

SNRK antibody (Biotin) (orb687263)

Description	SNRK antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human SNF-related serine/threonine-protein kinase protein (340-360AA)
Target	SNRK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NRH2
Expiration Date	12 months from date of receipt.