

Product Datasheet

IVNS1ABP antibody (Biotin) (orb687303)

Description	IVNS1ABP antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Influenza virus NS1A-binding protein (401-642AA)
Target	IVNS1ABP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6Y0
Expiration Date	12 months from date of receipt.