



Product Datasheet

ARL2BP antibody (orb687304)

Description	ARL2BP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 2-binding protein (1-163AA)
Target	ARL2BP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2Y0
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 20-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human glomerular tissue.



Western blot analysis of HepG2 whole cell lysate.