

Product Datasheet

plac8.2 antibody (Biotin) (orb687347)

Description	plac8.2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Zebrafish Zgc:158845 protein (5-19AA)
Target	plac8.2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A2VD52
Expiration Date	12 months from date of receipt.