

Product Datasheet

EXOC8 antibody (Biotin) (orb687423)

Description	EXOC8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Exocyst complex component 8 protein (508-725AA)
Target	EXOC8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYI6
Expiration Date	12 months from date of receipt.