

Product Datasheet

TMEM131L antibody (Biotin) (orb687435)

Description	TMEM131L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Transmembrane protein 131-like protein (963-978AA)
Target	TMEM131L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A2VDJ0
Expiration Date	12 months from date of receipt.