

Product Datasheet

CARD16 antibody (Biotin) (orb687571)

Description	CARD16 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Caspase recruitment domain-containing protein 16 protein (1-197AA)
Target	CARD16
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5EG05
Expiration Date	12 months from date of receipt.