


Product Datasheet

DNER antibody (orb687588)

Description	DNER antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Peptide sequence from Human Delta and Notch-like epidermal growth factor-related receptor protein (150-168AA)
Target	DNER
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NFT8
Dilution Range	IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U87 cells...