

Product Datasheet

CFHR5 antibody (FITC) (orb687682)

Description	CFHR5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Complement factor H-related protein 5 protein (374-569AA)
Target	CFHR5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXR6
Expiration Date	12 months from date of receipt.