

Product Datasheet

APOC2 antibody (Biotin) (orb687711)

Description	APOC2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Apolipoprotein C-II protein (20-38AA)
Target	APOC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02655
Expiration Date	12 months from date of receipt.