

Product Datasheet

TJP1 antibody (FITC) (orb687782)

Description	TJP1 antibody (FITC)
Species/Host	Other
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Tight junction protein ZO-1 protein (116-414aa)
Target	TJP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07157
Expiration Date	12 months from date of receipt.