


Product Datasheet

ANGPTL4 antibody (orb687808)

Description	ANGPTL4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Trefoil factor 3 protein (26-406AA)
Target	ANGPTL4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BY76
Dilution Range	WB: 1: 500-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y whole c...