



Product Datasheet

RPIA antibody (orb687820)

Description	RPIA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ribose-5-phosphate isomerase protein (1-311AA)
Target	RPIA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49247
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 200-1: 500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Western blot analysis of HepG2 whole cel...