


Product Datasheet

APOM antibody (orb687912)

Description	APOM antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Apolipoprotein M protein (1-188aa)
Target	APOM
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95445
Dilution Range	IHC-P: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of human
ki...