

Product Datasheet

CLCN6 antibody (Biotin) (orb687971)

Description	CLCN6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Chloride transport protein 6 protein (572-869AA)
Target	CLCN6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51797
Expiration Date	12 months from date of receipt.