


Product Datasheet

GLB1L3 antibody (orb688024)

Description	GLB1L3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Beta-galactosidase-1-like protein 3 protein (1-236AA)
Target	GLB1L3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NCI6
Dilution Range	WB: 1: 1000-1: 5000
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 whole cell...