



Product Datasheet

NFKBIA antibody (orb688108)


Description	NFKBIA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human NF-kappa-B inhibitor alpha protein (27-317AA)
Target	NFKBIA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25963
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 100-1: 500, IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...



Immunohistochemical staining of human ad...



Immunofluorescence analysis of HepG2 cel...