

Product Datasheet

OR4D10 antibody (HRP) (orb688165)

Description	OR4D10 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Olfactory receptor 4D10 protein (86-98AA)
Target	OR4D10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGI6
Expiration Date	12 months from date of receipt.