

Product Datasheet

OR5K2 antibody (HRP) (orb688177)

Description	OR5K2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Olfactory receptor 5K2 protein (6-24AA)
Target	OR5K2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NHB8
Expiration Date	12 months from date of receipt.