

Product Datasheet

OR5R1 antibody (Biotin) (orb688183)

Description	OR5R1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human Olfactory receptor 5R1 protein (292-307AA)
Target	OR5R1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NH85
Expiration Date	12 months from date of receipt.