

Product Datasheet

OR5T2 antibody (FITC) (orb688186)

Description	OR5T2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Olfactory receptor 5T2 protein (298-310AA)
Target	OR5T2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NKG2
Expiration Date	12 months from date of receipt.