



Product Datasheet

OR6B3 antibody (orb688192)


Description	OR6B3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Peptide sequence from Human Olfactory receptor 6B3 protein (259-272AA)
Target	OR6B3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGW1
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of U251 whole cell lysate.



Immunohistochemical staining of human skin tissue.