

Product Datasheet

OR6B3 antibody (FITC) (orb688194)

Description	OR6B3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Olfactory receptor 6B3 protein (259-272AA)
Target	OR6B3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGW1
Expiration Date	12 months from date of receipt.