


Product Datasheet

OSCAR antibody (orb688200)

Description	OSCAR antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence from Human Osteoclast-associated immunoglobulin-like receptor protein (265-279AA)
Target	OSCAR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYS5
Dilution Range	WB: 1: 1000-1: 5000
Expiration Date	12 months from date of receipt.



Western blot analysis of THP-1 whole cell...