


Product Datasheet

PCDHA10 antibody (orb688212)

Description	PCDHA10 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Protocadherin alpha-10 protein (276-324AA)
Target	PCDHA10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5I2
Dilution Range	IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y C...