

Product Datasheet

SDHB antibody (FITC) (orb688278)

Description	SDHB antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial protein (29-280AA)
Target	SDHB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21912
Expiration Date	12 months from date of receipt.