

Product Datasheet

TAS2R5 antibody (FITC) (orb688310)

Description	TAS2R5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Taste receptor type 2 member 5 protein (149-164AA)
Target	TAS2R5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYW4
Expiration Date	12 months from date of receipt.