


Product Datasheet

TBC1D3G antibody (orb688316)

Description	TBC1D3G antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human TBC1 domain family member 3G protein (312-549AA)
Target	TBC1D3G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6DHY5
Dilution Range	IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...