

Product Datasheet

TBC1D3G antibody (Biotin) (orb688323)

Description	TBC1D3G antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human TBC1 domain family member 3G protein (312-549AA)
Target	TBC1D3G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6DHY5
Expiration Date	12 months from date of receipt.