


Product Datasheet

GAPDH antibody (orb688347)

Description	GAPDH antibody
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep, Yeast
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic Peptide
Target	GAPDH
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04406
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa (1), Rat b...