


Product Datasheet

GAPDH antibody (orb688358)

Description	GAPDH antibody
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	GAPDH
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04406
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela using GAPD...